

AGISEAL<sup>®</sup>

# ELECTROSURGICAL GENERATOR

INTELLIGENT GENERATOR + REAL-TIME FEEDBACK TECHNOLOGY

## SL110M



(For oversea use only!)

SL110M




Security and Reliable



Intelligent Feedback



Convenient Operation



[ Vessel Sealing ]

Four Key Features

- Real-Time Monitoring, Safe Operation
- Intelligent Regulation, Efficient Hemostasis
- Low-Temperature Protection, Favorable Prognosis
- Multi-Purpose Instrument, Convenient and Reliable

Generator Parameters

| Weight           | Depth            | Width             | Height          |
|------------------|------------------|-------------------|-----------------|
| 5 kg (11.02 lbs) | 40 cm (15.75 in) | 32.5 cm (12.8 in) | 13 cm (5.12 in) |

Operating Output Characteristics

| Mode    | Rated Load | Peak Voltage | Rated Output Power | Operating Frequency |
|---------|------------|--------------|--------------------|---------------------|
| Sealing | 20Ω        | 287.5V       | 250W               | 472kHz ±10%         |

Nanjing ShouLiang-med Technology Co.,Ltd.

Tel:025-58370102 E-mail:info@shouliang-med.com

Block D9, No. 9 Kechuang Avenue, Intelligent Manufacturing Industrial Park (Zhihe Park),  
Jiangbei New Area, 211500 Nanjing City, Jiangsu Province, PEOPLE'S REPUBLIC OF CHINA

www.shouliangmed.com

