

Vessel Sealer Divider
2 major advantages


SEALING WITH FASTER SPEED

Used to safely seal large blood vessels, lymphatic vessels, and tissue bundles containing blood vessels or lymphatic vessels to form a permanent, clear, elastic sealing band.



Sealing speed as fast as **1-4 seconds**
1-3 times faster compared to traditional energy devices

(* According to our internal test report)



LOW TEMPERATURE JAW PROFILE
FASTER COOLING

Lateral thermal damage < 2mm
Less charring, sticking and smoke

(* According to our internal test report)

[Specifications Overview]

Model	Handle Type	Shaft Type	Shaft Length	Shaft Diameter	Shaft Rotation
Open curved big jaw series					
SL1220	Locking Handle	Rectangular Shaft	18cm	13.8mm	180°
SL1320EA	Press-type Handle	Rectangular Shaft	18cm	13.8mm	360°
SL1420EA	Press-type Handle	Round Shaft	18cm	14.8mm	360°

Model	Handle Type	Jaw Curvature	Overall Length	Seal Length	Cut Length
Open scissors series					
SL5300S	Scissors Press	28°	19cm	16.5mm	14.7mm
SL1500EA	Scissors Press	40°	21cm	20.6mm	19.8mm
SL1600EA	Scissors Press	40°	25cm	24.5mm	21.8mm
SL1700EA	Scissors Press	28°	22.8cm	23.9mm	21.4mm



[Applicable Surgeries]

Including but not limited to the following departments



Head and Neck Surgery

Thyroidectomy, Parotidectomy,
Neck lymph node dissection



General Surgery

Gastrectomy, Colorectal resection,
Pancreatectomy, Mastectomy



Gynecology

Open hysterectomy, Vaginal hysterectomy,
Adnexectomy



Urology

Nephrectomy, Cystectomy,
Prostatectomy

Nanjing ShouLiang-med Technology Co.,Ltd.

Tel: 025-58370102 E-mail: info@shouliang-med.com
D9 block, NO 9 Kechuang avenue, Jiangbei new area, Intelligent Manufacturing
Industrial Park (Zhihe Park), Nanjing city, Jiangsu, China.

SHOULIANG-MED

www.shouliangmed.com

ELECTROSURGICAL
VESSEL SEALER / DIVIDER


OPEN SURGERY
COOLER SEALING BLOODLESS CUTTING

(For overseas use only !)

Brand-new open curved big jaw series

Unique long seal length for efficient dissection; newly designed curved handle for ergonomic, lightweight operation

Model:SL1320EA



360°continuous rotation for flexible access to any surgical field.



Press-type energy activation, sealing and cutting in 2 steps, simplifying operation.



Soft-touch gel grip with four anti-slip zones for comfortable long-term use.

Other models:



SL1420EA



SL1220

Rectangular Shaft

Large diameter shaft for stable control

Double-sided Thermal Insulation

90°constant temperature with multi-layer protection to reduce thermal injury

Advanced Design Features

Nano-coating to reduce tissue adhesion; textured surface for enhanced grip;ceramic dots for insulation and short-circuit prevention



14°Curved Jaw

Maryland curved jaws fit organ contours; 36 mm atraumatic seal length ensures precise, efficient open surgery.



Brand-new open scissors series

Multiple configurations for all open surgical scenarios; Stepwise tactile feedback system to optimize precision.



(For overseas use only !)

SHOULIANG-MED

